

DIG-labeling of RNA Probes for In Situ Hybridization

by Tatsuya Okafuji, Kevin Mitchell Lab, Smurfit Institute of Genetics, Trinity College Dublin

Updated on 20/10/2009

1. Add the following reagents in an RNase-free tube. Mix gently by pipetting.

10x transcription buffer	2 μ l
Template DNA	1 μ g (plasmid) or 0.25 μ g (PCR product)
10x DIG labeling mix (Roche)	2 μ l
RNase inhibitor (40 U/ μ l)	0.5 μ l
T7, T3 or SP6 RNA polymerase (*1)	20U
Nuclease-free H ₂ O (Sigma)	to 20 μ l

2. Incubate at 37°C for 2-4 hr. Prepare an agarose gel for electrophoresis.

3. After 1.5 or 3.5 h*, take a 1 μ l aliquot to run on an agarose gel for quality check of RNA.

4. Add 1 μ l of RNase-free DNase I and incubate for 15 min at 37°C.

5. Stop reaction by adding 2 μ l of 0.2 M EDTA pH 8.0.

6. Column-purify RNA [We use NucleoSpin RNA Clean-up XS (Macherey-Nagel) and elute RNA in 30 μ l of RNase-free H₂O.

7. Quantify RNA (We usually do that with Nanodrop using 1.5 μ l of eluted RNA).

8. Adjust the concentration of RNA to 0.1 μ g/ μ l with DEPC-treated H₂O.

9. Add 1/100 volume of 100 mM sodium citrate pH 6.35.

10. Store at -20°C until use.

* RNA probes can be checked after finishing all the procedures.